

# NEW JETQUICK PCR Purification Protocol

## VACUUM MANIFOLD PROTOCOL INSERT

September 2003

### Very Important!

Before starting the procedure, make sure that solution H2 is reconstituted as indicated on the bottle's label.

Removal of mineral oil is **not** required.

The JETQUICK procedure is **not** sensitive to detergents (i.e. Triton X-100, Nonidet NP40) or gelatin in the assay. Nonetheless, final concentrations of >0.1% of Tween-20 should be avoided.

All centrifugation steps are carried out at top speed (approx. 13,000 rpm) in a conventional table-top microfuge.

### PROTOCOL

- 1.) **Sample Preparation:** Add **400 µl** of binding **buffer H1** to **100 µl of amplification reaction** and **mix thoroughly**.

***Important note:** Make sure, that the volume of the PCR assay does not exceed 100 µl (WITHOUT any oil overlay). If PCR assays with a volume >100 µl are to be processed, remove most of the oil overlay with a pipette, and adjust the volume of binding buffer H1 proportionally to a ratio of buffer H1 : amplification reaction of 4:1 (v/v). In such a case multiple successive loading steps (see step 3) are necessary.*

- 2.) **Vacuum Manifold Preparation:** Attach the vacuum manifold to a vacuum source. Attach as many micro-spin columns as required to the female luer extensions on the vacuum manifold.
- 3.) **Column Loading:** Load the sample from step 1 into a prepared micro-spin column. Apply vacuum (-200 to -650 mbar) until all liquid has been pulled through the cartridge. Then turn off the vacuum source.
- 4.) **Column Wash:** Add **700 µl of reconstituted Wash Buffer H2** to the micro-spin column. Apply vacuum (-600 to -650 mbar) until all liquid has been pulled through the cartridge. Then turn off the vacuum source.  
Place the micro-spin column into a 2 ml wash tube and centrifuge at top speed for 1 min to remove any residual wash buffer.
- 5.) **DNA Elution:** Place the JETQUICK micro-spin column into a new 1.5 ml recovery tube and add **50 µl** of sterile water (or TE buffer or 10 mM Tris-HCl [pH 7-8]) directly onto the center of the silica matrix of the spin column. Centrifuge at top speed for 2 min.

***Important:** Higher DNA concentrations will be obtained if the elution is carried out with only 30 µl of elution buffer. In this case preheat the elution buffer to 65-70 °C, add the buffer onto the center of the silica matrix of the spin column and let stand for 1 min before centrifugation. DNA eluted in water should be generally stored at -20 °C.*